## **XAS** CABINETS up to $60^{\circ C}$



The XAS range is used to dry and store glassware by injecting air into the working volume at an adjustable temperature up to 60°C. They are functional, easy-to-use and -to-clean, and energy efficient.



www.france-etuves.com

Fixed feet.

Height adjustable.

# XAS RANGE GLASSWARE DRYING CABINETS

Swivel castors with brakes.

Height adjustable. Easy movement.

#### **TECHNICAL DATA**



	XAS320	XAS800	
SIZES			
Working volume (L)	320	800	
Internal width (mm)	900	900	
Internal depth (mm)	600	600	
Internal height (mm)	600	1480	
External width (mm)	946	946	
External depth (mm)	713	741	
External height (mm)	878	1836	
Weight, empty (kg)	92	115	
Number of door(s)	2	2	
PERFORMANCES			
Temperature range (°C)	up to 60°C		
ELECTRICAL DATA			

ELEC	CTRICAL DATA		
Heating power (W)	2000	2000	
Power	230V ~1P 50/60Hz, 16A		
LOADING SHELVES			
Maximum number	6	12	
Maximum distributed load per shelf (kg)	30	30	
Maximum load in the cabinet (kg)	60	120	

#### **ACCESSORIES**



Shelf with guide bars

- Shelf with guide bars
  - Stainless steel. Perforated. Maximum distributed load: 30 kg.
- 2 Subframe with casters (for XAS320) Elevates the cabinet for a easier loading. Swivel castors with brakes. Shelf surcharge.

Subframe with castors





#### Order a few click...

#### www.france-etuves-store.com

Web, mobile, tablet, order easily and receive your products quickly!



### www.france-etuves-store.com





65/67 rue Auguste Meunier 77500 Chelles, FRANCE

T: +33 1 64 21 30 60 / F: +33 1 64 21 30 86