

# Ultrasonic Pipet Washer

## AW-31

Pipet Qty.

1mL  
264pcs

5mL  
216pcs

10mL  
136pcs

High Frequency  
Output

55W

- With return current system for smart cleaning efficiency.
- Dirty liquid is drained by the force of jet from bottom nozzle for good cleaning.
- Bath basket is spinning constantly for good cleaning.
- Compact design.
- Washing and rinsing can be done in one washer.
- Washable by water-soluble alkaline detergent.

### Specifications

Product code	291704
<b>Model</b>	<b>AW-31</b>
Pipet quantity	Pipet (Max. 500mm)
	1mL = 264pcs.
	5mL = 216pcs.
	10mL = 136pcs.
Operation water amount	1 to 3L/min.
Ultrasonic oscillator	High frequency output, 55W, 28kHz
	Oscillation plate: stainless steel, SUS430
	Oscillator: ball ferrite
Bath material	Clear vinyl chloride
Cabinet	Stainless Steel SUS304
Rack dimensions (Basket)	I.D. $\phi$ 138×H540mm
	O.D. $\phi$ 146×H570mm
External dimensions*1	W390×D400×H844mm
Weight	Main body: approx. 18.5kg, Basket: approx. 1.0kg
Power (50/60Hz)	AC115V/AC220V Single phase with step-down transformer
Accessories	Rack
	Manual pump
	Water supply hose : I.D. 12mm×1.5m length
	Drain hose : I.D. 25mm×1.5m length
	Detergent

\*1 Do not include protrusions.

